

Aseptic Containment Isolator (ACTI)



Aseptic Containment Isolator (ACTI)

Introduction

- Esco Aseptic Containment Isolator (ACTI) work in conjunction with advanced material transfer techniques and bio decontamination agents providing a 6 log reduction in viable contaminants.
- Esco ACTIs provide standard configurable designs able to adapt to various batch sizes and process flows. Through a fully user selectable operating system, the same Isolator can be setup to operate under recirculation or full exhaust airflow and operate in positive or negative pressure modes, allowing the system to be multifunctional and cater for all requirements of toxic or non-toxic aseptic materials. For toxic materials, the system incorporates safe change filters.


Basic principles

- Full unidirectional airflow provides superior aseptic work zones.
- Safety toughened laminated glass hinges upwards assisted with gas springs for batch loading.
- Airflow regime runs either in recirculating or full exhaust airflow. Full exhaust airflow for fast purging of bio decontamination agent during aeration period and recirculation option for reduced airflow taken from the room and exhausted during normal operation and during conditioning and decontamination phases.
- U15 supply and exhaust filters suitable for either toxic or non-toxic aseptic materials. Exhaust Filter directly below the isolator can be removed & a bypass tube fitted for non-potent material applications.

Standard Features

- Full unidirectional airflow provides superior aseptic work zones.
- Safety toughened laminated glass hinges upwards assisted with gas springs for batch loading.

Airflow Regimes



GENERAL SPECIFICATIONS

		ACTI-2G_	ACTI-3G_	ACTI-4G_
Nominal Size Main		1.2m	1.6m	2.0m
Working Chamber Dimensions		1.2m x 0.75m x 0.85	1.6m x 0.75m x 0.85	2.0m x 0.75m x 0.85
External Dimensions (W x D x H)	With Adjustable Base Stand (Min)	1.2m x 1.1m x 2.2m	1.6m x 1.1m x 2.2m	2.0m x 1.1m x 2.2m
	With Adjustable Base Stand (Max)	1.2m x 0.75m x 2.68m	1.6m x 0.75m x 2.68m	2.0m x 0.75m x 2.68m
Glove Port Height Min		1080	1080	1080
Glove Port Height Max		1360	1360	1360
Chamber Environment		ISO Class 5 all Chambers (Grade A)		
Filter Type – Isolator Inlet		ULPA U15 with Integral mesh guard and knife edge gel seal		
Filter Efficiency - Inlet		99.9998%	99.9998%	99.9998%
Filter Type – Isolator Exhaust		HEPA H14 with knife edge gel seal		
Filter Efficiency - Exhaust		99.997%	99.997%	99.997%
Lighting Level		≥600Lux	≥600Lux	≥600Lux
Sound Level		≤68dBA	≤68dBA	≤68dBA
Isolator Construction	Chamber	316L	316L	316L
	Service Housing	316L	316L	316L
	Support Frame	316L	316L	316L
Isolator Finish	Chamber Internal	≤0.4Ra	≤0.4Ra	≤0.4Ra
	Chamber External	≤0.6Ra	≤0.6Ra	≤0.6Ra
	Service Housing External	0.6Ra	0.6Ra	0.6Ra
	Support Frame	1.0Ra	1.0Ra	1.0Ra
Electrical Requirements (By Client)	220-240V, AC, 50Hz, 1Ø	✓	✓	✓
	110-120V, AC, 60Hz, 1Ø	✓	✓	✓
	220-240V, AC, 60Hz, 1Ø	✓	✓	✓
Compressed Air requirements (By Client)	6 Barg Pressure at 5ltr/sec	✓	✓	✓
Exhaust Duct requirements (By Client)		10" Duct from Isolator to Outside		
Options	Pass Chamber	✓	✓	✓
	Bio-Decontamination Steris	✓	✓	✓
	Bio-Decontamination Bioquell	✓	✓	✓
	Non-Viable Air Sampler	✓	✓	✓
	Viable Air Sampler	✓	✓	✓
	Sterility Test Pump	✓	✓	✓
	Glove Tester	✓	✓	✓
	Waste Bag Grommet	✓	✓	✓
	Sterile Continuous Liner	✓	✓	✓
	Bag-Welder with Table	✓	✓	✓
	RTP Ø270 - Alpha	✓	✓	✓
	RTP Ø270 - Beta Canister	✓	✓	✓
	RTP Ø270 - Beta Liner	✓	✓	✓

ESCO PHARMA
PLATFORM SPECIALIST.

Esco Pharma Pte Ltd
21 Changi South Street 1 Singapore 486777 • Tel: +65 65420833
Fax: +65 65426920 • Email: csis.pharma@escoglobal.com

Esco Technologies, Inc.
1661 Loretta Ave., Feasterville, PA 19053
Tel: 215 322 2155 • Email: eti.pharma@escoglobal.com

Esco Gb Ltd
Unit 9, Century Park Networkcentre, Dearne Lane,
Manvers, Rotherham, South Yorkshire, S63 5DE.
Tel: +44 (0) 1709 761 669 • Email: egb.info@escoglobal.com

escoglobal.com
LIVE SUPPORT

ISOCIDE™

Esco Global Offices: Bahrain | Bangladesh | China | India | Indonesia | Italy | Japan
Malaysia | Philippines | Russia | Singapore | South Africa | South Korea | Thailand
United Kingdom | USA Vietnam



Esco Pharma ACTI, Solihull, A4, Feb/2016
 Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials. Esco
 reserves the right to alter its products and specifications without notice. All trademarks and logos in this material
 are the property of Esco and the respective companies.

PT Esco Britain Indonesia
Cert. No: 851333E

Esco Micro Pte Ltd
Cert. No: 851276

PT Esco Britain Indonesia
Cert. No: QN0508070001