

Streamline[®] Compounding Isolator

Compounding Aseptic Containment Isolators
(Recirculating and Total Exhaust/Single Pass)



SCI-2G8-N1RL-1-0-EG



SCI-3G8-N1SL-1-0-SS




STREAMLINE®
COMPOUNDING ISOLATOR

Streamline® Compounding Isolator (SCI) provides an aseptic environment for pharmaceutical preparations. It can be configured to operate in a recirculating or total exhaust air flow scheme.


Also, aside from positive pressure, this particular unit has two negative pressurization classifications, the purposes of which are:

- N1: factory-configured to provide a more negative pressured pass-through chamber (PTC) in accordance with the CETA guideline for testing compounding isolators.
- N2: this unit is intended for powder handling applications, thus, the work chamber (minimum negative -37 Pa) is more negative than the PTC (minimum negative -25 Pa).

Overall, the SCI is designed to provide user, product, environment, and cross-contamination protection from exposure to hazards.



SCI Stainless Steel Recirculating Unit



SCI Electrogalvanized Steel Recirculating Unit

Applications

Pharmacy Compounding (Hazardous / Non-hazardous)	Potent Formulations	Sterility Testing (Small Scale)	TPN Formulation	Hazardous Drugs Compounding
	Contained Powder Handling	Cell Culture	Gene and Cell Therapy	Tissue Engineering

Main Features

- Class 3 Leak Tight Containment, as per ISO 10648-2.



SCI Electrogalvanized Steel Total Exhaust Unit

Safety and Certification

	Design	Cabinet Performance	Air Cleanliness	Electrical Safety
Standard Compliance	USP <797> and <800>, USA NIOSH, OSHA	Class 3 Leak Tight Containment as per ISO 10648-2, CETA CAG-002-2006	ISO 14644-1, Class 5, EU GMP Grade A, Worldwide JIS B9920, Class 3, Japan BS 5295, Class 1, UK	IEC 61010-1, Worldwide EN 61010-1, USA CAN/ CSA-22.2, No. 61010-1



SCI Total Exhaust Unit



Main Features

- ULPA filters (as per IEST-RP-CC001.3 and HEPA (H14) filter as per EN 1822) with a typical efficiency of >99.995% at Most Penetrating Particle Size (MPPS) 0.1-0.3 μm , providing ISO Class 5 air cleanliness as per ISO 14644-1.
- For SCI with Exhaust Type 1 and 3, Bag-In Bag-Out (BIBO) Exhaust Filter at the top allows for safe and convenient filter change.
- Sentinel™ Gold Microprocessor Controller supervises all functions and monitors airflow and pressure in real-time.
- For negative pressure configuration, the work zone and passthrough interchange are under negative pressure to the room to maintain operator protection in the event of a breach in the isolation system.
- For positive pressure configuration, both the work zone and pass-through chambers are under positive pressure to the room to ensure product protection.
- Robust dual-wall construction. Unique Esco Dynamic Chamber™ plenum surrounds filter seals with negative pressure.
- Ergonomically angled front to improve reach and operator comfort.
- Adaptor plate for utilities such as: Sharps Bin Beside the Work Zone, Service Fixtures, Electrical Outlets, etc.
- Safe-change cuff rings permit glove changes with zero risk of contaminating the work zone or environment.
 - Sleeves are made of Hypalon material while gloves are determined by user requirements.
- The unit has FDA-grade air-tight static/compression seals.
- It has easy to clean work zone without any crevices.
- Options for external material of construction:
 - Electrogalvanized steel with ISOCIDE™ antimicrobial powder coating
 - Full stainless steel 304
- A hydraulic height-adjustable stand is available as an option.
- The unit is equipped with Foot Switch to open the inner doors.
- It has options to have a sharps disposal system:
 - bottom of the work zone
 - beside the work zone




SCI Recirculating Unit





The United States Pharmacopeia (USP) General Chapters <797> and <800> together with the FDA 503B Guidance, states the need for proper documentation of the standard operating procedures (SOPs) surrounding the activities for sterile pharmaceutical preparations.

- Utilizes energy-efficient light-emitting diode (LED) which has higher durability and better lighting quality.
- Internal chamber of SCI consists of two (2) planes made from single-piece SS 316L with 17.5 mm radius coved corners.
- Work top with a drain pan at the bottom.
- Electromagnetic interlocking door with timedelayed ingress/egress control to minimize particle entry; assuring work zone sterility during material transfer.
- Optional connectivity for external exhaust fan with zero volt relay interlock.



Esco's integration options for a **CCTV** and a **monitor** and **keyboard mount** can easily provide real-time monitoring and recording of personnel activities; meeting the demands of international guidelines.

CCTV integration:


- front panel with a stainless steel mount
- rear wall of the work zone

Monitor and Keyboard integration:


- access to rear view monitor system
- beside the front panel equipped with a keyboard and mouse



Adjustable foot rest



IV bar with hooks



CCTV



Audio-Visual Alarm Package
for Sentinel Microprocessor
Controller




Granite Slab

- SS 316 Frame
- Leveling Feet

AIRFLOW PATTERNS FOR STREAMLINE® COMPOUNDING ISOLATOR MODELS

Recirculating Airflow

- Ambient air is pulled through the inlet pre-filter G4 and supply filter located on top of the isolator.
- The downflow filter creates a laminar airflow providing ISO Class 5 air cleanliness (ISO 14644-1) to the main chamber and pass-through, with a down flow velocity of 0.40 m/s +/- 20%.
- Air from the work zone and pass-through is quickly purged out by the exhaust fan to keep the area clean. The fan pulls approximately 70% of the purged air back to the plenum, passes through the HEPA downflow filter again, and is recirculated back to the work zone and pass-through.
- Approximately 30% of the purged air is exhausted through the filter to prevent heat build-up inside the isolator.
- Exhausted air is replenished by ambient air coming from the top inlet pre-filter G4 that has 80% efficiency.



- HEPA-filtered air
- Unfiltered / potentially contaminated air
- Room air / Inflow air

Guide to Streamline® Compounding Isolator Models


SCI-3G8-N1R

Isolator Unit	Model	No. of gloves - Nominal Width		Voltage Code		Pressure	
Streamline® Compounding Isolator	SCI	2G	5 ft (1.6m)	8	220-240 V AC, 50/60Hz, 1Ø	P0	Process zone: + ≥ 37 Pa (min) Pass-through zone: + ≥ 25 Pa (min)
		3G	6 ft (1.95m)	9	110-120 V AC, 50/60Hz, 1Ø	N1	Process zone: - 25 Pa (min) Pass-through zone: - 37 Pa (min)
						N2	Process zone: - 37 Pa (min) Pass-through zone: - 25 Pa (min)

Upon ordering, input material of construction at the end of the model code: SCI-2GC_-N1RL-1-SE-SS or -EG

- SS: Full stainless steel exterior
- EG: Electrogalvanized steel with ISOCIDE™ coating

* The configuration allows for 1 sharps container and 1 stainless steel vessel for other laboratory waste.



- HEPA-filtered air
- Unfiltered / potentially contaminated air
- Room air / Inflow air

Total Exhaust/Single Pass

- Ambient air is pulled through the inlet pre-filter G4 and supply filter via the main fans at the top of the isolator. This creates positive pressure on the plenum which provides the downflow of air with a velocity of 0.40 m/s \pm 20%.
- The pre-filter extends the life of the filters by trapping larger particulates that can easily clog the main filters.
- This downflow supply then provides an ISO Class 5 environment and unidirectional airflow inside the isolator; thus, protecting the materials inside the main chamber and pass-through.
- Air from the work zone and pass-through is then quickly purged by the fans to keep the area clean. The purge is completely exhausted through HEPA filters as well; ensuring that only clean air is exhausted back to the environment.

ounding Isolator (SCI) models

N1RL-1-SE-SS

Airflow Scheme		Pass-Through zone		Exhaust Type	Code	Sharps Provision		External Construction	
R	Recirculating	L	1 Left	Double HEPA Exhaust (BIBO Top and non-BIBO Bottom)	1	Without	0	Electrogalvanized steel	EG
S	Total Exhaust/ Single Pass	R	1 Right	Single Exhaust (Top, non-BIBO)	2	With (below work zone) 5 gal: 190.50 x 266.70 x 431.80mm (7.5" x 10.50" x 17.0")	SB	Stainless steel	SS
				Single Exhaust (Top, BIBO)	3	With (beside work zone) 5 gal: 190.50 x 266.70 x 431.80mm (7.5" x 10.50" x 17.0")*	SE		

GENERAL SPECIFICATIONS
Streamline® Compounding Isolator


GENERAL SPECIFICATIONS Streamline® Compounding Isolator			SCI-2G _ _	SCI-3G _ _	
Unit Nominal Size (Width)			1645 mm (5')	1950 mm (6')	
External Dimension (W x D x H) - with 1 left pass chamber*	Exhaust Type - 2 (Single Exhaust, Top, Non-BIBO)	Without stand	1645 x 825 x 1330 mm (64.76 x 33.27 x 52.36")	1950 x 825 x 1330 mm (76.77 x 33.27 x 52.36")	
		With SPC-A/ SPL-A (720mm, 28")	1645 x 825 x 2050 mm (64.76 x 33.27 x 80.71")	1950 x 825 x 2050 mm (76.77 x 33.27 x 80.71")	
		With SPC-B/ SPL-B (860mm, 34")	1645 x 825 x 2190 mm (64.76 x 33.27 x 86.22")	1950 x 825 x 2190 mm (76.77 x 33.27 x 86.22")	
		With STL (660 to 960 mm)	1645 x 825 x 1990 to 2290 mm (64.76 x 33.27 x 78.35 to 90.16")	1950 x 825 x 1990 to 2290 mm (76.77 x 33.27 x 78.35 to 90.16")	
		With SHM (685 to 935 mm)	1645 x 825 x 2015 to 2265 mm (64.76 x 33.27 x 79.33 to 89.20")	1950 x 825 x 2015 to 2265 mm (76.77 x 33.27 x 79.33 to 89.20")	
	Exhaust Type - 1 (Double Exhaust, Top BIBO, Below Non-BIBO)**	Without stand	1645 x 825 x 1650 mm (64.76 x 33.27 x 64.57")	1950 x 825 x 1650 (76.77 x 33.27 x 64.57")	
		With SPC-A/ SPL-A (720mm, 28")	1645 x 825 x 2370 mm (64.76 x 33.27 x 92.91")	1950 x 825 x 2370 mm (76.77 x 33.27 x 92.91")	
		With SPC-B/ SPL-B (860mm, 34")	1645 x 825 x 2510 mm (64.76 x 33.27 x 98.43")	1950 x 825 x 2510 mm (76.77 x 33.27 x 98.43")	
		With STL (660 to 960 mm)	1645 x 825 x 2310 to 2610 mm (64.76 x 33.27 x 90.55 to 102.36")	1950 x 825 x 2310 to 2610 mm (76.77 x 33.27 x 90.55 to 102.36")	
		With SHM (685 to 935 mm)	1645 x 825 x 2335 to 2580 mm (64.76 x 33.27 x 91.54 to 101.18")	1950 x 825 x 2335 to 2580 mm (76.77 x 33.27 x 91.54 to 101.18")	
Process Chamber Internal Dimension (W x D x H)			960 x 625 x 625 mm (37.80 x 24.61 x 24.61")	1265 x 625 x 625 mm (49.80 x 24.61 x 24.61")	
Pass Chamber Internal Dimension (W x D x H)			450 x 625 x 625 mm (17.72 x 24.61 x 24.61")		
Pass Chamber Opening Dimension (W x H)		Inner Door	315 x 445 mm (12.40 x 17.52")		
		Outer Door	355 x 445 mm (13.98 x 17.52")		
Pass Chamber and Process Chamber Door Material			Polycarbonate		
Airflow Regime			Factory Configured Recirculating or Total Exhaust/Single Pass		
Pressurization			Factory Configured Positive or Negative Pressure		
Glove Port Diameter			200 mm (Circular)		
Glove Port Quantity			2	3	
Chamber Environment			ISO Class 5 for all chambers (Grade A)		
Process Chamber Downflow Velocity			0.4 +/- 20% m/s (1.31 fps)		
Pre-filter			G4, panel, polyester fiber media		
Downflow and Exhaust Filter Type			HEPA (H14) Filter with Integral Metal Guards and Filter Frame Gaskets; Fully Compliant With EN 1822 (H14) and IEST-RPCC001.3 Requirements		
Filter Efficiency			>99.995% for particle 0.1-0.3 microns (MPPS, as per EN1822)		
Lighting Level			>650 Lux		
Sound Level			≤ 67 dBA		
Isolator Construction	Main Body		1.2 mm (0.05") 18 gauge electro-galvanized steel with white oven-baked epoxy-polyester antimicrobial powdercoated finish*		
	Work Tray		1.5 mm (0.06") 16 gauge stainless steel, type 316L, with 4B finish		
	Side Walls		1.2 mm (0.05") 18 gauge stainless steel, type 316L, with 4B finish		
Electrical	220-240V, AC, 50 Hz, 1Ø		SCI-2G8-_-_-_-_-_-_-_-	SCI-3G8-_-_-_-_-_-_-_-	
	Cabinet Full Load Amps (FLA)		10 A	11 A	
	Optional Outlets FLA		5 A	5 A	
	Cabinet Nominal Power		238 W, 1.5 A	333 W, 2.3 A	
	Cabinet BTU		812 BTU/hr	1136 BTU/hr	
	110-120V, AC, 50 Hz, 1Ø		SCI-2G9-_-_-_-_-_-_-_-	SCI-3G9-_-_-_-_-_-_-_-	
	Cabinet Full Load Amps (FLA)		Contact Esco for more details		
	Optional Outlets FLA				
	Cabinet Nominal Power				
Cabinet BTU					
Options	Carbon Filter		Contact Esco for more details		
	Circular glove port		5180346		
	CCTV Provision		5180046		
	Drain		Contact Esco for more details		
	Electrical Outlet		Contact Esco for more details		
	Alarm package		5170227		
	Pre-filter, G4		5090114		
	Glove Leak Tester - Circular		With client-supplied compressed air: 5180311 / With Mobile Compressor: 5180312		
	IV Bar with SS hooks		5170930	5170931	
	UV Lamp		5170251	5170255	
	Rear View Screen Adaptation		5180033		
	Anti-blow Back Valve (ABBV)		EG-Steel: 5170352 / Stainless Steel: 5170354		
	Additional Manual Latches		5180038		
	Top Exhaust Collar		EG-Steel: 5171251 / Stainless Steel: 5171253		
	Side Exhaust Collar		EG-Steel: 5171252 / Stainless Steel: 5171254		
	Automated Pressure Hold Test		With client-supplied compressed air: 5180031 / With Mobile Compressor: 5180032		
	Sharps Disposal Container, 18.7 L		Positive pressure: 5170223 / Negative pressure: 5170224		
	Shipping Weight			500 kg (1102.31 lbs)	700 kg (1543.68 lbs)
	Shipping Dimensions, Maximum (W x D x H)			1720 x 900 x 2250 mm (67.72 x 35.43 x 88.58")	2100 x 1150 x 2250 mm (82.68 x 45.28 x 88.58")

Note:
*All common surgical gloves can be attached to the cuff rings; thereby, making the system adaptable to operator preference.
*No different of height for the unit with and without sharps bin disposal below the workzone
*To add ABBV, please add +420mm for the overall height
*For Exhaust Type-2, to add Exhaust Carbon Filter, please add +65 mm for the overall height



*For Exhaust Type-2, to add Exhaust Collar, please add +105mm for the overall height
*For Exhaust Type-1 and 3, to add Top Exhaust collar with/without Exhaust Carbon Filter, please add +205mm for the overall height
*For Exhaust Type-1 and 3, to add Side Exhaust collar with/without Exhaust Carbon Filter, please add +165mm for the overall height
** For Exhaust Type-1 and 3, the standard base unit will need to be connected to ducting.

ENGINEERING DRAWING

SCI Recirculating Model (SCI-2G8-N1RL-2-0-SS)



19



1

2


3

4

5

6

1645



7

8

9

10

11

12


1330

2050

720

625*

825



13

14

15






16

17

*625

- | | |
|--|---|
| 1. Esco Sentinel Microprocessor Control System | 11. UV Tube 15w |
| 2. Circular Glove Ports 200x200mm | 12. Pass Chamber Polycarbonate Inner Door |
| 3. Main Chamber Polycarbonate Door | 13. Exhaust Fan |
| 4. Support Stand SPC-5A0-S-SS G2 | 14. Pass Chamber Supply Fan |
| 5. Pass Chamber Polycarbonate Outer Door | 15. Pass Chamber Supply HEPA H14 Filter |
| 6. Inner Door Foot Switch | 16. Electrical Outlet Code C |
| 7. Exhaust HEPA H14 Filter | 17. Pass Chamber Sliding Tray |
| 8. Main Chamber Supply Fan | 18. LED Light |
| 9. Main Chamber Supply HEPA H14 Filter | 19. G4 Air Inlet Pre-Filter |
| 10. IV Bar With Hooks (6pcs) | |





ENGINEERING DRAWING
SCI Recirculating Model (SCI-3G8-P0RL-2-SE-SS)



- | | |
|--|---|
| 1. Esco Sentinel Microprocessor Control System | 11. IV Bar With Hooks (6pcs) |
| 2. Circular Glove Ports 200x200mm | 12. UV Tube 30w |
| 3. Main Chamber Glass Door | 13. Pass Chamber Glass Inner Door |
| 4. Sharp Disposal Container (Red) | 14. LED Light |
| 5. Motorized Hyd S/Stand Sh | 15. Exhaust Fan |
| 6. Pass Chamber Glass Outer Door | 16. Pass Chamber Supply Fan |
| 7. Inner Door Foot Switch | 17. Pass Chamber Supply HEPA H14 Filter |
| 8. Exhaust HEPA H14 Filter | 18. Electrical Outlet Code C |
| 9. Main Chamber Supply Fan | 19. Pass Chamber Sliding Tray |
| 10. Main Chamber Supply HEPA H14 Filter | 20. G4 Air Inlet Pre-Filter |

ENGINEERING DRAWING

SCI Total Exhaust Model (SCI-2G8-N2SL-1-0-EG)



- | | |
|--|---|
| 1. Esco Sentinel Microprocessor Control System | 12. Main Chamber Supply HEPA H14 Filter |
| 2. Circular Glove Ports 200x200mm | 13. IV Bar With Hooks 6pcs |
| 3. Main Chamber Polycarbonate Door | 14. UV Tube 15w |
| 4. Support Stand SPC-5A0-S-G2 | 15. Pass Chamber Polycarbonate Inner Door |
| 5. Exhaust Collar (Top of unit) | 16. LED Light |
| 6. Pass Chamber Polycarbonate Outer Door | 17. Pass Chamber Supply Fan |
| 7. Inner Door Foot Switch | 18. Pass Chamber Supply HEPA H14 Filter |
| 8. Exhaust Carbon Filter | 19. Electrical Outlet (Provision Only) |
| 9. Exhaust HEPA Filter, Bag-In Bag-Out (BIBO) | 20. Pass Chamber Sliding Tray |
| 10. Exhaust Fan | 21. G4-Air Inlet Pre-Filter |
| 11. Main Chamber Supply Fan | |

Isolation Technology

Isolation containment systems provide inherently superior sterility compared to open front clean air devices such as laminar flow clean benches and Class II biological safety cabinets. USP <797> guidelines specify that isolators may be situated in an area subject to less severe environmental controls compared with open front clean air devices.

Isolators are an effective solution especially for lower-volume pharmaceutical processes. They reduce operating and renovation costs, take up less space, and are easier to maintain. Positive pressure isolators are designed for non-hazardous aseptic processes, while negative pressure isolators are intended to promote personnel and environment safety during processing of hazardous materials such as antineoplastic or cytotoxic compounds.


- For volatile hazardous drugs, it is recommended to have a negatively pressured, total exhaust isolator. Alternatively, a recirculating isolator with optional exhaust activated carbon filter and/or external ventilation, can be selected.

Maximum Patient Protection and Product Sterility

The combination of downflow and exhaust HEPA filters provides the Esco Streamline® Compounding Isolator, a fully integrated envelope for product and operator protection in all configurations.

- Advanced mini-pleated supply HEPA filters are tested to >99.995% efficiency for to 0.3 micron particulates; significantly better than conventional filters.
 - An improved mini-pleat separation technique maximizes filter surface area, improves efficiency and extends filter life over conventional separation.
 - The HEPA supply filter provides an ISO Class 5 air in the work surface via a gentle vertical laminar flow.
 - All filters have sampling ports for filter integrity testing during qualification.
- Laminar airflow within the work zone and pass-through enables recovery of chamber atmosphere to ISO Class 5 conditions within 3 minutes following a worse-case contamination event. The entire work zone air is changed 20-30 times per minute.

Mini-pleat Separatorless Filter (left) vs. Conventional Aluminum Separator Filter (right)




Esco cabinets use Swedish Camfil Farr® mini-pleat filters without aluminum separators to increase filter efficiency, minimize the chance of leakage, and prolong filter life. Filters include a lightweight aluminum frame for structural stability and elimination of swelling common to conventional wood frames.

- Optional sharps disposal system enables a smoother work flow while minimizing material transfer to enhance patient protection and product sterility.
 - The interface between sharps disposal bin and isolator is aerosol tight to avoid ingress of contamination during the disposal operation.
- Safe-change cuff rings enable glove change with zero risk of contamination.
- The exterior construction of an isolator can either be of full stainless steel or electrogalvanized steel with an ISOCIDE™ antimicrobial powder coating to minimize contamination.
 - Performance results are available upon request.

Esco Filter Efficiency

(%) Typical Penetration



● Typical Penetration

The SCI unit uses ULPA filters (as per IEST-RP-CC001.3 and HEPA (H14) filter as per EN 1822) with a typical efficiency of >99.995% at 0.1 to 0.3 microns, providing ISO Class 5 air cleanliness as per ISO 14644-1.

Fan Efficiency

The Streamline® Compounding Isolator fan system is designed for maximum energy efficiency with minimal maintenance.

- Designed with a centrifugal, direct-drive, and external rotor motors to reduce operating costs.
- Unique Esco motor/fan orientations minimize noise and vibration.
- Built-in solid-state variable speed controllers are infinitely adjustable from Off to Maximum




Electrical Outlets and Utility Fittings

- Electrical outlet, ground fault, North America
- UV lamp
- IV bar with hooks
- Latches
- Damper
- Alarm package
- Exhaust Collar
- Duplex and/or GFI outlet
- Perforated shelf to increase work zone space
- Sharps bin container
- Side exhaust connection (for Exhaust Type 1 and 3)
- Back-up battery for electromagnetic interlocking doors (Contact Esco for more information)
- Granite Slab
- Anti Blowback Valve

Highlights

- Output 10 VDC, max. 1.1 mA, Tach output, Control input 0-10 VDC, PWM
- Motor current limit, soft start, locked-rotor protection, over temperature protected electronics/motor
- Condensate discharge holes: none, open rotor
- Ingress protection of IP54
 - Material: Fiberglass reinforced composite
 - Impeller: PA6
 - Direction of rotation: Clockwise, seen on rotor

ISOCIDE™ Antimicrobial Powder-Coating



All exterior painted surfaces are powder-coated with Esco Isocide™, an antimicrobial inhibitor to minimize contamination. Isocide™ is integrated into the coating substrate and cannot wash out or diminish by repeated cleaning. Performance results are available upon request. Contact Esco or your Esco Sales Representative for details.

Sentinel™ Microprocessor Control, Monitoring System

The Esco Sentinel™ microprocessor-based control system supervises operation of all cabinet functions. Controls are configurable to meet user requirements.


Enlarged, multi-line digital read-out with alpha-numeric display indicates all input, status and alarm functions.

Work zone, interchange pressures and downflow velocity are continuously measured and displayed. Integrated true airflow sensor is temperature compensated for improved accuracy.

Password-protected administration can be set to restrict access to main menu.

Color coded indicator lamps display green for fan operation; blue for LED lights and electrical outlet; and orange for AUX function ON.

Optional audible and visual alarms warn of unsafe conditions such as low chamber pressure.



Sentinel Microprocessor Control System, Programmable

- When programmed ON • the start-up sequence confirms status with Air Safe and local time display.
- the Personal Identification Number (PIN) access restricts unauthorized adjustments.

Main Control Panel Features

- Work zone and pass-through pressures are monitored and displayed on the LCD screen
- Continuous monitoring and display of cabinet laminar airflow on large, easy-to-read LCD display.
- An optional alarm package (pressure and airflow) is available for users with more sophisticated requirements.

Safe Glove Change Procedure: Replacing Disposable Gloves

Safe change design system allows glove change at the middle of a process or when the equipment is in operation.



1. Pull the Glove/Sleeve outside the isolator.



2. Fold the fingers of the glove inside the cuff ring.



3. Remove the outer ring.



4. Carefully roll the gloves from the middle groove to the outer groove.



5. Take the new glove and ensure the thumb is at the top. Stretch the ring of the new glove over the port and over the old glove onto the middle groove.



6. Install the ring up to the middle groove.



7. Carefully loosen the old glove from the outer groove.



8. Put the glove/sleeve inside the isolator.



9. Working with one hand in the adjacent glove, carefully pull the old glove.



10. The procedure is now complete.

Safe Sleeve Change Procedure: Replacing the Sleeves



1. Remove the screws that secure the glove port cover



2. Remove the outer glove port cover



3. Remove the "O" ring



4. Carefully roll the ring of the sleeves/gloves from the inner groove to the outer groove of the port



5. Ensure that the old sleeves/gloves is inside the isolator



6. Take the new sleeves and ensure the thumb is at the top and stretch the "O" ring of the new sleeves over the port and over the old sleeves into the inner groove



7. Replace the "O" ring into the outer groove of the glove port



8. Working with one hand in the adjacent sleeves, carefully work from the outer ring and into the isolator. The old sleeves needs to be remove while under the new sleeves




9. Return the glove port outer cover.



10. Secure the port cover with the screws. The procedure is now complete

ESCO LIFESCIENCES GROUP NETWORK

42 Locations in 21 Countries All Over the World



- Global Offices
- Distributors
- Factories
- R&D Centers
- Regional Distribution Centers



Air Shower
 Aseptic Containment Isolator (ACTI)
 Ceiling Laminar Airflow Units
 Cleanroom Transfer Hatch
 Containment Barrier Isolator (CBI)
 Downflow Booth (DFB)
 Dynamic Floor Laminar Hatch
 Dynamic Pass Box
 Evidence Drying Cabinet
 Garment Storage Cabinet
 General Processing Platform Isolator (GPPI)
 Laminar Flow Horizontal Trolley
 Laminar Flow Straddle Units, Single and Double
 Laminar Flow Vertical Trolley
 Pass Box
 Soft Wall Cleanroom
 Sputum Booth
 Ventilated Balance Enclosure (VBE)
 Weighing and Dispensing Containment Isolator (WDCI)

Since 1978, Esco has emerged as a leader in the development of controlled environment, laboratory and pharmaceutical equipment solutions. Products sold in more than 100 countries include biological safety cabinets, fume hoods, ductless fume hoods, laminar flow clean benches, animal containment workstations, cytotoxic cabinets, hospital pharmacy isolators, and PCR cabinets and instrumentation. With the most extensive product line in the industry, Esco has passed more tests, in more languages, for more certifications, throughout more countries than any biosafety cabinet manufacturer in the world. Esco remains dedicated to delivering innovative solutions for the clinical, life science, research and industrial laboratory community.

ESCO
 HEALTHCARE



Esco Healthcare
 19 Changi South Street 1,
 Singapore 486779
 Tel: +65 65420833
 Email: mail@vaccixcell.com

Esco Technologies, Inc.
 2512 Metropolitan Drive, Suite 120 B
 Feasterville- Trevose, PA 19053-6738
 Tel: +1 215 322 2155
 Email: eti.pharma@escolifesciences.com

Esco GB Ltd
 Unit 2 R-evolution @ Gateway 36, Kestrel
 Way, Barnsley, S70 5SZ
 Tel: +44 (0) 1226 360 799
 Email: egb.info@escolifesciences.co

Esco Lifesciences Offices: Bangladesh | China | Denmark | Germany | Hong Kong | Indonesia | Italy | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam



9010607_Streamline@ Compounding Isolator Brochure_A4_vA_Blue_03042025
 Esco can accept no responsibility for possible errors in catalogues, brochures and other printed materials.
 Esco reserves the right to alter its products and specifications without notice. All trademarks and logotypes in this material are the property of Esco and the respective companies.