

CLAF

Spare Parts



Consumable Parts Listing

Tag Number	Part Name	Description	Order Code
Downflow Filter	Filter, H14	Megalam MG 915 x 730 x 105 mm	1100381
XL01 to XL06	LED Tear Drop Light 3 ft	LED T5 Lamp, 86-256 VAC, 16 W, 3 ft., 3000 k, 50/60 Hz	N/A

Accessories and Spare Parts Listing

Tag Number	Part Name	Description	Order Code
MAE01 to MAE04	Blower ECM Centrifugal	380-480 V/50/60 Hz,3,2A, D3G-250-EE51-11	1110053
PDT1101	Pressure Differential Transmitter	4-20 mA, 0-500 Pa	1281621
PDT1102	Pressure Differential Transmitter	4-20 mA, -250/+250 Pa	1281622
PS001	Power Supply Block	24 VDC Power Supply, 180 - 550 VAC, 24 VDC @ 5A WR - 120 - 24	1281635
SI03	Light Indicator	Light indicator, Red, 230 VAC	1281688
SI02	Light Indicator	Techna LED TEC, Yellow, 230 VAC	1281707
EA1	Light Indicator	Schneider, LED, 24 VDC, Red	1080763
EA2	Light Indicator	Schneider, LED, 24 VDC, Green	1080762
KF1	DC relay for lightnings	24 VDC coil, RXM3AB2BD	1281612
KF2 to KF6	DC relay for alarm dry-contacts	24 VDC coil, PLC-BSC-24DC/21-21	1281517
FLT01-FLT02	Outlet filter for control panel	IP 54 including filter mat for ventilation by convection	1281515