

DFB

Spare Parts



Consumable Parts Listing

Part Name	Description	Order Code
Downflow Filter H14	Megalam MG 715 x 730 x 105	1100379
Pre-filter	Aeropleat 30/30 - 594 x 594 x 44 mm	1100362
F8 Filter	591 mm x 591 mm x 292 mm	1100361
T5 Lamp	FLFT5EH24W865	1081004

Accessories and Spare Parts Listing

Part Name	Description	Order Code
EC Fan	3Ph 380-480 v 50/60 Hz, 355 mm DIA	1281617
Pressure Differential Transmitter	4-20 mA, 0-750 Pa, CX8MB2421KPA	1282067
Pressure Differential Transmitter	4-20 mA, 0-500 Pa, CX8MB24500PA	1282066
Temperature and Relative Humidity Transmitter	4-20 mA, RPHP3E44LCD	1282068
Electronic Ballast	WH5-230-L	1081003