

VBE

Spare Parts



Consumable Parts Listing

Part Name	Description	Order Code
HEPA Filter	457 x 762 x 90	1100029
Bag, PVC change out	80C x 96-3G w/ 16" clear	1287017
O-ring, 234 x 232 x 6 mm	Thickness, 40 shore	1287111

Spare Parts Listing

Part Name	Description	Order Code
Relay Board	Board, 5 Relay, C-Series	1080025
Main Board	Board, Main, C-Series	1090539
Power Supply	SMPS Meanwell PT65D	1080305
Membrane	Membrane Grey	1080247
LCD	Display, LCD, 1 Line	1080154
Interface Board	Board, Interface, C-Series	1080027
Sensor	Transmitter, Pressure	1281621
Lamp Holder	Holder, Lamp, T5	1080205
Fluorescent Tube	Lamp, Fluorescent, 24 W	1081004
Ballast	Fluorescent Ballast 120 V, 0.56 A, 50/60 Hz	1080009
Buzzer	Buzzer, 400 mm lead	1080038
Resistor	Wirewound, 500 ohm, 1%	1281609
Speed Control	RS485 to 0-10 V Convertor	1081130
Circuit Breaker	Circuit Breaker, 6A	1080119
Blower	Blower, EC, Centrifugal (BLR-38)	1110028