

WDCI

Spare Parts



Category	Tag Number	Description	Range	Connection	Quantity	Esco Code Number
Filters	HF101, HF102, HF103	Push-Push Filter	35m ³ /hr at 130 Pa at 99.98% efficiency	Push-Push	3	1100325
	HF201, HF202	Filter for Glove Boxes - Closed Casing PVC	30m ³ /hr at 180 Pa at 99.98% efficiency	Ø61.6, Thread 2.5	2	1100326
	HF104	Instrument Filter	Acro 50 Vent Device with PTFE Membrane	Hose Barb	2	1050771
	PF101, PF102	Pre-Filter - Synthetic Media	1.1 m/s at 41 Pa at 90% efficiency	Ø150 mm Disc, 1" thick	4	1100327
Valves	BFVP101, BFVP102	Butterfly Valve - Pneumatic - NO	Pneumatically actuated, indicator, and limit switch	3" Triclover - BS4825-3	1	1288002
	BFVP201, BFVP202					
Gloves/ Glove Ports	GL 101 to GL 105 GP 201	Gauntlets	10", Hypalon, 15/24 Thickness, 32" long, 9.5" hand size	O' Ring	1	1286005
	GP 101 to GP 105 GP 201	Glove port cover	200 mm x 300 mm elliptical	Push Fit	2	1286002
	GP 101 to GP 105 GP 201	Glove port ring	200 mm x 300 mm elliptical	Bolted	2	1286003
	GP 101 to GP 105 GP 201	Glove port backing ring	200 mm x 300 mm elliptical	Bolted	2	1286004
Gaskets and Seals		WDCI 2G Inflatable Seal - Visor Door	2190 profile (16 mm x 16 mm) 7 mm inflation	Inflatable Seal	1	1287147
		WDCI 3G Inflatable Seal - Visor Door	SP - ESC - 009	Inflatable Seal	1	1287148
		Outer Gasket - Glove port	2 mm thick	Bolted	2	1287003
		Inner Gasket - Glove port	3 mm thick	Bolted	2	1287004
		Inflatable Seal Drum Head Size #1	SP - ESC - 006	Inflatable Seal	1	1287149
		Inflatable Seal Drum Head Size #2	SP - ESC - 007	Inflatable Seal	1	1287150
		Inflatable Seal Drum Head Size #3	SP - ESC - 008	Inflatable Seal	1	1287151
		WDCI 2G RTP cover O-ring	OD 104 ID 93 x 6 mm thickness 40 shore	Static Seal	1	1288404
		WDCI 3G RTP cover O-ring	OD 234 x 222 mm 6 mm thickness 40 shore	Static Seal	1	1288405

Category	Tag Number	Description	Range	Connection	Quantity	Esco Code Number
Mechanical		Gas Struts - Variable Pressure	250 mm stroke x 545 mm overall	M8 x 1.25 Thread	1	1288042
		Castor	Ø5" wheel	Plate	4	1288031
		1"/1.5" Triclover Clamp	N/A	Single Hinge - Wingnut	4	1288021
		1" Triclover Seal	N/A	N/A	4	1288020
		1" Triclover Blank	N/A	N/A	2	1288022
		3" Triclover Clamp	N/A	Single Hinge - Wingnut	4	1288024
		3" Triclover Seal	N/A	N/A	4	1288025
		4" Triclover Clamp	N/A	Single Hinge - Wingnut	2	1288026
		4" Triclover Seal	N/A	N/A	2	1288027
		4" Triclover Blank	N/A	N/A	1	1288028
Fans and Motors	FA101	Centrifugal Fan c/w Direct Coupled Motor	Inverter Ready, motor fitted with Thermoresistors - 0.37 Kw	N/A	1	1289003
	FA 201					
Electrical	XL101	Fluorescent Light with enclosure	220-240 VAC, 50 Hz, T5 2 x 24 W	120-277 VAC, 50 to 60 Hz	1	1281418
	XL201				1	
	ZES101	Proximity Switch	4 mm operating distance	8 mm Sensing Range	1	1281403
	ZES201				1	
	ZES 102	Distance Measurement Sensor	4-20 mA, 5-30 VDC supply	0.2 to 8 m	1	1281401
	INT 101	Electromagnetic lock	120 lbs	2 Wire, 12 VDC	1	1080157
	INT 201				1	
		LS101, LS102, LS201, LS202	Dual Sensor (Normally Open)	3 mm sensing distance, 6-60 VDC	6-60 VDC	1
Instrumentation	PDT101	Differential Pressure Transmitter	-300 Pa to +300 Pa	Horse Barb - For 6 mm ID Tube	1	1281415
	PDT201	Differential Pressure Transmitter	0 to 500 Pa	Horse Barb - For 6 mm ID Tube	1	1281417
	PDT102	Differential Pressure Transmitter	0 to 500 Pa	Horse Barb - For 6 mm ID Tube	1	1281417
	PDT103	Differential Pressure Transmitter	0 to 500 Pa	Horse Barb - For 6 mm ID Tube	1	1281417
	PDT104	Differential Pressure Transmitter	0 to 500 Pa	Horse Barb - For 6 mm ID Tube	1	1281417
	PDT105	Differential Pressure Transmitter	0 to 500 Pa	Horse Barb - For 6 mm ID Tube	1	1281417
	PDT202, PDT2013	Differential Pressure Transmitter	0 to 300 Pa	Horse Barb - For 6 mm ID Tube	1	1281416
Pneumatics	PS101	Pressure Switch	15-30 VDC, 0-10 bar		1	1288017
	PS102 PS201 PS202	Pressure Switch	15-30 VDC, 0-2 bar		1	1288011
	REG0101	Filter regulator	3/8" BSP		1	1288014